

according to Regulation (EC) No 1907/2006

Lab Putty Base

Revision date: 20.04.2015

Print date: 30.04.2015

Page 1 of 3

SECTION 1: Identification of the substance/mixture and of the company/undertaking**1.1. Product identifier**

Lab Putty Base

1.2. Relevant identified uses of the substance or mixture and uses advised against**Use of the substance/mixture**

for dental use only

1.3. Details of the supplier of the safety data sheet

Company name:	Coltène/Whaledent AG
Street:	Feldwiesenstrasse 20
Place:	CH-9450 Altstätten
Telephone:	+41 (71) 75 75 300
Telefax:	+41 (71) 75 75 301
e-mail:	info.ch@coltene.com
Internet:	www.coltene.com
Responsible Department:	Regulatory Affairs msds@coltene.com

Further Information

Only supplied to dentists and dental laboratories or upon their instructions

SECTION 2: Hazards identification**2.1. Classification of the substance or mixture**

This mixture is classified as not hazardous according to Regulation (EC) 1272/2008 [GHS].
The preparation is not dangerous in the sense of Directive 1999/45/EC.
Preparation/mixture without hazards for human health or environment

2.2. Label elements**Hazardous components which must be listed on the label**

According to EC directive 1999/45/EC or annex VI of directive 67/548/EEC the preparation does not have to be labelled.

SECTION 3: Composition/information on ingredients**3.2. Mixtures****Chemical characterization**

Polymer preparations and compounds

SECTION 4: First aid measures**4.1. Description of first aid measures****General information**

No special measures are necessary.

After contact with eyes

Rinse immediately carefully and thoroughly with eye-bath or water.

SECTION 5: Firefighting measures**5.1. Extinguishing media****Suitable extinguishing media**Water. Foam. Extinguishing powder. Carbon dioxide (CO₂). Sand.

according to Regulation (EC) No 1907/2006

Lab Putty Base

Revision date: 20.04.2015

Print date: 30.04.2015

Page 2 of 3

Additional information

Co-ordinate fire-fighting measures to the fire surroundings.

SECTION 6: Accidental release measures**6.2. Environmental precautions**

Do not allow to enter into surface water or drains.

6.3. Methods and material for containment and cleaning up

Wipe up with absorbent material (eg. cloth, fleece). Clear contaminated areas thoroughly.

SECTION 7: Handling and storage**7.1. Precautions for safe handling****Advice on safe handling**No special precautionary measures are necessary.
Observe instructions for use.**7.2. Conditions for safe storage, including any incompatibilities****Further information on storage conditions**Keep the packing dry and well sealed to prevent contamination and absorption of humidity.
Recommended storage temperature: 15 - 23 °C**SECTION 8: Exposure controls/personal protection****8.1. Control parameters****8.2. Exposure controls****Protective and hygiene measures**

When using do not eat, drink or smoke.

Hand protectionSingle-use gloves.
Unsuitable material: NR (Natural rubber (Caoutchouc), Natural latex).**SECTION 9: Physical and chemical properties****9.1. Information on basic physical and chemical properties**

Physical state: Paste

Colour:

Test method**Changes in the physical state**

Initial boiling point and boiling range: >150 °C

Water solubility: insoluble
(at 23 °C)**Solubility in other solvents**

ketone.

SECTION 10: Stability and reactivity**Further information**

No special measures are necessary.

SECTION 11: Toxicological information**11.1. Information on toxicological effects**

according to Regulation (EC) No 1907/2006

Lab Putty Base

Revision date: 20.04.2015

Print date: 30.04.2015

Page 3 of 3

Acute toxicity

LD50: >2000 mg/kg (Rat)

The statement is derived from products of similar composition.

Irritation and corrosivity

Irritant effect on the skin: (Rabbit)

Not an irritant.

The statement is derived from products of similar composition.

Sensitising effects

Guinea-pig.

no danger of sensitization.

The statement is derived from products of similar composition.

SECTION 12: Ecological information**Further information**

No special environmental measures are necessary.

SECTION 13: Disposal considerations**13.1. Waste treatment methods****Advice on disposal**

Can be burnt together with household waste in compliance with official regulations in contact with approved waste disposal companies and with authorities in charge.

SECTION 14: Transport information**Other applicable information**

Not a hazardous material with respect to these transportation regulations.

SECTION 15: Regulatory information**15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture****National regulatory information****Additional information**

SECTION 16: Other information**Further Information**

The above information describes exclusively the safety requirements of the product and is based on our present-day knowledge. The information is intended to give you advice about the safe handling of the product named in this safety data sheet, for storage, processing, transport and disposal.

(The data for the hazardous ingredients were taken respectively from the last version of the sub-contractor's safety data sheet.)